



Fluid Path Products

Safety. Repeatability. American Made.



Gauge Guards



Why FPP Gauge Guards?

Gauge Guards protect your critical instrumentation from process chemicals, and can substantially extend their lifespan. Our high-performing polymer and elastomer compounds provide longer run times and less contamination than competing brands.

Sizes and Materials

Gauge guards are available in sizes 1/2"-3" and in materials such as EPDM, Platinum Cured Silicone, and Teflon.

Our Materials Meet the Following Requirements:

- FDA Extraction Studies per 21CFR177.2600
- Testing according to USP Class VI, Chapter <87> and Chapter <88>
- ADI Free
- RoHS Directive 2011/65/EU and 2015/863

Fluid Path Products customizes! Our in-house tooling design can typically take your ideas from concept to finished product within two to three weeks at a cheaper cost than our competitors.

Size	Platinum Silicone	EPDM	FKM	PTFE
3/4"	42GG-PMSE-075	42GG-EZ-075	42GG-SFY-075	42GG-TN-075
1"	42GG-PMSE-100	42GG-EZ-100	42GG-SFY-100	42GG-TN-100
1 1/2"	40GG-PMSE-150	40GG-EZ-150	40GG-SFY-150	40GG-TN-150
2"	40GG-PMSE-200	40GG-EZ-200	40GG-SFY-200	40GG-TN-200
3"	40GG-PMSE-300	40GG-EZ-300	40GG-SFY-300	40GG-TN-300

